

- O Reduces risk of accidental displacement of catheter
- O Safe anchoring and protection of catheter
- O Prevents contamination of puncture site from environment
- O Allows inspection of puncture site
- Reduces catheter movement
- O Can prevent kinking of catheter

Drain-Guard 8150

Anchoring device for percutaneous drainage catheters

Drain-Guard

Safe anchoring of catheters

Drain-Guard is an anchoring device for percutaneous drainage catheters that are protruding from the patient's skin in an angle approximately perpendicular to the patient's skin surface. The device consists of a skin-side surface with a skin friendly adhesive combined with a specially designed foam block that safely anchors catheters from 5 to 16 French.


The soft foam support block can prevent kinking of the catheter. A transparent hydrocolloid cover allows continuous monitoring of the puncture site and provides an efficient protection against contamination.

For the sake of the patient

Drain-Guard's soft and pliable adhesive carrier can offer comfort to the patient. Because the hydrocolloid has been carefully tested for skin irritation, the risk of allergic reactions is minimal. The Drain-Guard gives the patient great freedom to move around without worrying about the cathe-

Easy to use

The Drain-Guard catheter anchoring device is easy to apply to the patient's skin and changing of the device may even be carried out by a family member.



Ordering information:

Drain-Guard



REF.	Product	Packing
8150	Drain-Guard - Achoring device for percutaneous drainage catheters.	1 unit in peel-pouch, sterile 25 units per box

Distributor:



EN ISO 13485:2016 Quality Management System

©2021 Copenhagen MedLab ApS. All rights reserved.

Copenhagen MedLab ApS

- Formervangen 5 DK-2600 Glostrup Danemark
- Tel.: +45 43 29 66 88 Fax: +45 43 29 66 77
- info@copenhagen-medlab.com